

Omnium1 – Bundles

Omnium1 iMRS one Basic Set

Omnium1 Kit
OmniMat
OmniPad
D/A Converter
iMRS one Android software



Omnium1 iMRS one Complete Set

Omnium1 Kit
OmniMat
OmniPad
OmniSpot
D/A Converter
iMRS one Android software



Omnium1 iMRS one Basic Combo

Omnium1 Kit
OmniMat
OmniPad
OmniBrain
D/A Converter
iMRS one android software



Omnium1 iMRS one Complete Combo

Omnium1 Kit
OmniMat
OmniPad
OmniSpot
OmniBrain
D/A Converter
iMRS one Android software



Omnium 1 OmniBrain Set

Omnium1 Kit
OmniBrain-Hardware
OmniBrain Android software




What does „PEMF“ stand for?

The term PEMF means „pulsed electromagnetic fields“ and describes in modern medicine the specific application of unique pulsed, electromagnetic frequencies and intensities for the health benefit of human beings and animals.

Research has shown that low energy pulsating fields within the frequency range of “the biological window” (3-25 Hz), as well as the inner and outer wave patterns of the Earth’s magnetic field (7.83-11.79 Hz) can produce remarkable positive changes within the body. In exploring this technology for your health and the health of your family, insist upon true PEMF with the iMRS one W@pp!

The term “magnetic resonance stimulation” simply indicates that by applying a pulsating electromagnetic field that matches the Earth’s magnetic field, nearly all 75 trillion cells of the body can be stimulated into harmonic resonance. As a result the powerful capacity within the cell membrane increases and enables the cells to produce more energy. This energy enhances the self-healing powers of the body.



250-833-8781

www.joyfulhealthstore.omnium1.com

www.omnium1.com

© Swiss Bionic Solutions Holding GmbH

OMNIUM1
INTELLIGENT LIFESTYLE



OWN YOURSELF!

Omnium1 – Health goes Purple

Developed, engineered and designed in Switzerland the Omnium1 is based on the worldwide known and established Android platform. It can be used as a regular tablet without any restrictions: Surfing the web via WIFI, sending, receiving emails and complete access to the android app structure. The brand new android version 4.4.2 KitKat is already preinstalled and serves as the user interface for all related applications. In contrast to a regular market-based tablet the Omnium1 not only contains of a specific high performance battery (14,8V, 9000mA/h), it is to date the world`s first tablet equipped with two serial ports (specific 20-pin connector). These accesses allow connectivity and control of a wide range of health, wellness and fitness related hard- and software applications (W@pp`s). Daily routines like blood pressure evaluation, blood glucose measuring, ECG`s, EEG`s, heart rate variability and even active and analog applications like pulsed electro magnetic field therapy (PEMF) or brain wave entrainment can soon be performed, evaluated and monitored by and with only ONE single control unit.



Prevention, active health, wellness and fitness management in the comfort of your home is no longer a vision. Due to its versatility and performance the Omnium1 provides comprehensive home health management in an easy, slick and cost effective way.

www.pemf-joy.com

iMRS one

www.myjoyfulhealth.com

The iMRS one-W@pp consists of a software package, an intelligent D/A converter to transform the digital control signals into analog, pulsating, electromagnetic fields (PEMF) and a hardware-set (applicators) for the output (multiple options available). iMRS one is the world`s first magnetic resonance stimulations system W@pp which performs on an android platform. It is technically based on the most sold and most popular PEMF system in the world for home use, the iMRS.

Together with the Omnium1 interface the iMRS one consists of the unique capability of performing socket free, easy to use and totally independent PEMF applications with various applicators. The software surface is very intuitive, easy to operate and due to the proper size of the touch icons also suitable for the elderly. The new generation of applicators (OmniMat, OmniPad and OmniSpot) are outstanding, due to their high quality fabric, excellent workmanship and a compact construction for an easy and smooth transportation.

OmniBrain

Over decades Brain Wave Entrainment has become an integral part of modern mental training concepts. The technology is widely used in modern medicine as well as in professional sports and as a very effective, supportive learning tool. The primary goal of a brain wave entrainment session is to utilize the basic tendency of the human brain to resonate with external stimuli. Hereby one can apply visible or audible, external frequencies which can be recognized and processed by the brain. By controlling the applied frequencies in a proper way, the brain starts resonating with it and adjusts its own frequency patterns. This way, the user can actively influence the actual mental state. Omnium1 is a 3-way-brain wave entrainment system, which is capable of gaining resonances based on light, tone and color. The goggle is made of suede-leather, very comfortable to wear and contains of built in high end ear phones for the tone and music stimulus as well as a flexible PCB with RGB-LED`s to create the desired light and color frequencies.

